

## New Shamir Relax™

### Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Height	Invisible Marking
Hard Resin	1.498	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	14mm	15 Ⓢ RLX 5
Blue Zero	1.498	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	14mm	BS5 Ⓢ RLX 5
Polarized	1.498	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	14mm	PZ Ⓢ RLX 5
Transitions® Signature®	1.5	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray, Brown & Green	14mm	TR Ⓢ RLX 5
Transitions® XTRActive®	1.5	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray & Brown	14mm	X5 Ⓢ RLX 5
Transitions® Drivewear®	1.5	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Drivewear®	14mm	DW Ⓢ RLX 5
DLC™ (Trivex)	1.53	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	14mm	TX Ⓢ RLX 5
DLC™ (Trivex) Blue Zero	1.53	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	14mm	BZX Ⓢ RLX 5
DLC™ (Trivex) Polarized	1.53	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	14mm	PTX Ⓢ RLX 5
DLC™ (Trivex) Transitions® Signature®	1.53	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray, Brown & Green	14mm	TTX Ⓢ RLX 5
DLC™ (Trivex) Transitions® XTRActive®	1.53	-14.00 to +8.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Gray & Brown	14mm	XX Ⓢ RLX 5
1.56	1.56	-14.00 to +8.50	0.5, 2, 3, 4, 5, 6, 7, 8, 9	0.75 to 3.50	Clear	14mm	56 Ⓢ RLX 5

\* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



## New Shamir Relax™

### Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Height	Invisible Marking
Polycarbonate	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	14mm	PC Ⓢ RLX 5
Polycarbonate Blue Zero	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	14mm	BZP Ⓢ RLX 5
Polycarbonate Polarized	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	14mm	PPZ Ⓢ RLX 5
Polycarbonate Transitions® Signature®	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray, Brown & Green	14mm	PT Ⓢ RLX 5
Polycarbonate Transitions® XTRActive®	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray & Brown	14mm	PX Ⓢ RLX 5
Polycarbonate Transitions® Drivewear®	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Drivewear®	14mm	PD Ⓢ RLX 5
SuperLite™ 1.6	1.6	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	14mm	16 Ⓢ RLX 5
SuperLite Blue Zero 1.6	1.6	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	14mm	BZ6 Ⓢ RLX 5
SuperLite 1.6 Polarized	1.6	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Clear	14mm	PZ6 Ⓢ RLX 5
SuperLite 1.6 Transitions® Signature®	1.6	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray, Brown & Green	14mm	T6 Ⓢ RLX 5
SuperLite 1.6 Transitions® XTRActive®	1.6	-14.00 to +8.75	0.5, 2, 3, 4, 5, 6, 7, 8	0.75 to 3.50	Gray & Brown	14mm	X6 Ⓢ RLX 5
SuperLite 1.67	1.67	-16.00 to +16.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9.5, 11, 12.5	0.75 to 3.50	Clear	14mm	67 Ⓢ RLX 5

\* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



## New Shamir Relax™

### Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Height	Invisible Marking
SuperLite Blue Zero 1.67	1.67	-16.00 to +16.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9.5, 11, 12.5	0.75 to 3.50	Clear	14mm	BZ67 RLX 5
SuperLite 1.67 Polarized	1.67	-16.00 to +16.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9.5, 11, 12.5	0.75 to 3.50	Clear	14mm	PZ67 RLX 5
SuperLite 1.67 Transitions® Signature®	1.67	-16.00 to +16.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9.5, 11, 12.5	0.75 to 3.50	Gray, Brown & Green	14mm	T67 RLX 5
SuperLite 1.67 Transitions® XTRActive®	1.67	-16.00 to +16.00	0.5, 2, 3, 4, 5, 6, 7, 8, 9.5, 11, 12.5	0.75 to 3.50	Gray & Brown	14mm	X67 RLX 5
SuperLite 1.74	1.74	-20.00 to +16.75	0.5, 1, 1.5, 2.5, 3.5, 4.75, 6, 7, 8, 9, 9.75, 10.5, 11.25, 12	0.75 to 3.50	Clear	14mm	74 RLX 5
SuperLite 1.74 Transitions® Signature®	1.74	-20.00 to +16.75	0.5, 1, 1.5, 2.5, 3.5, 4.75, 6, 7, 8, 9, 9.75, 10.5, 11.25, 12	0.75 to 3.50	Gray & Brown	14mm	T74 RLX 5

\* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.



**SHAMIR** ReCreating Perfect Vision®