





New Shamir Relax™




Minimum Fitting Height: 14mm





New Shamir Relax™



Invisible Marking RLX 5	
Relax +0.50	RLX 5
Relax +0.65	RLX 6
Relax +0.80	RLX 8

Ink Marking	
Relax +0.50	Relax 5
Relax +0.65	Relax 6
Relax +0.80	Relax 8

Material	Invisible Markings (XX)
HR 1.5	15
Transitions® 1.5	TR
XTRActive® 1.5	X5
Polarized 1.5	PZ
Blue Zero 1.5	BZ5
Drivewear® 1.5	DW
Trivex® 1.53	TX
Transitions® 1.53	TTX
XTRActive® 1.53	XX
Polarized 1.53	PTX
Blue Zero 1.53	BZX
Mid Index 1.56	56
Polycarbonate 1.59	PC
Polycarbonate 1.59 Transitions®	PT
Polycarbonate 1.59 XTRActive®	PX

Material	Invisible Markings (XX)
Polycarbonate 1.59 Polarized	PPZ
Blue Zero 1.59	BZP
Drivewear 1.59	PD
SuperLite™ 1.6	16
SuperLite™ 1.6 Transitions®	T6
SuperLite™ 1.6 XTRActive®	X6
SuperLite™ 1.6 Polarized	PZ6
SuperLite™ 1.6 Blue Zero	BZ6
SuperLite™ 1.67	67
SuperLite™ 1.67 Transitions®	T67
SuperLite™ 1.67 XTRActive®	X67
SuperLite™ 1.67 Polarized	PZ67
SuperLite™ 1.67 Blue Zero	BZ67
SuperLite™ 1.74	74
SuperLite™ 1.74 Transitions®	T74

Ink Marking: Blue
Invisible Marking: Red

Fitting Point
Prism Reference Point