

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Height	Invisible Marking
DLC™ (Trivex)	1.53	-5.00 to +4.00	6, 7, 8	1.00 to 3.50	Clear	19mm	TX S 250
DLC™ (Trivex) Transitions®	1.53	-5.00 to +4.00	6, 7, 8	1.00 to 3.50	Gray & Brown	19mm	TTX S 250
DLC™ (Trivex) Transitions® XTRActive™	1.53	-5.00 to +4.00	6, 7, 8	1.00 to 3.50	Gray	19mm	XX S 250
DLC™ (Trivex) Polarized	1.53	-5.00 to +4.00	6, 7, 8	1.00 to 3.50	Gray & Brown	19mm	PTX S 250
1.56	1.56	-5.50 to +5.25	6, 7, 8	1.00 to 3.50	Clear	19mm	56 S 250
Polycarbonate	1.586	-6.00 to +5.25	6, 7, 8	1.00 to 3.50	Clear	19mm	PC S 250
Polycarbonate Transitions®	1.586	-6.00 to +5.25	6, 7, 8	1.00 to 3.50	Gray & Brown	19mm	PT S 250
Polycarbonate Transitions® XTRActive™	1.586	-6.00 to +5.25	6, 7, 8	1.00 to 3.50	Gray	19mm	PX S 250
Polycarbonate Polarized	1.586	-6.00 to +5.25	6, 7, 8	1.00 to 3.50	Gray & Brown	19mm	PPZ S 250

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Add	Color	Fitting Height	Invisible Marking
Superlite™ 1.6	1.6	-6.00 to +5.25	6, 7, 8	1.00 to 3.50	Clear	19mm	16 S 250
Superlite 1.6 Transitions®	1.6	-6.00 to +5.25	6, 7, 8	1.00 to 3.50	Gray & Brown	19mm	T6 S 250
Superlite 1.6 Transitions® XTRActive™	1.6	-6.00 to +5.25	6, 7, 8	1.00 to 3.50	Gray	19mm	X6 S 250
Superlite 1.6 Polarized	1.6	-6.00 to +5.25	6, 7, 8	1.00 to 3.50	Gray & Brown	19mm	PZ6 S 250

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Color	Invisible Marking
DLC™ (Trivex)	1.53	-6.00 to +6.00	6, 7, 8	Clear	TX SV
DLC™ (Trivex) Transitions®	1.53	-6.00 to +6.00	6, 7, 8	Gray & Brown	TTX SV
DLC™ (Trivex) Transitions® XTRActive™	1.53	-6.00 to +6.00	6, 7, 8	Gray	XX SV
DLC™ (Trivex) Polarized	1.53	-6.00 to +6.00	6, 7, 8	Gray & Brown	PTX SV
Polycarbonate	1.586	-6.00 to +6.00	6, 7, 8	Clear	PC SV
Polycarbonate Transitions®	1.586	-6.00 to +6.00	6, 7, 8	Gray & Brown	PT SV
Polycarbonate Transitions® XTRActive™	1.586	-6.00 to +6.00	6, 7, 8	Gray	PX SV
Polycarbonate Polarized	1.586	-6.00 to +6.00	6, 7, 8	Gray & Brown	PPZ SV
Superlite™ 1.6	1.6	-6.00 to +6.00	6, 7, 8	Clear	16 SV
Superlite 1.6 Transitions®	1.6	-6.00 to +6.00	6, 7, 8	Gray & Brown	T6 SV
Superlite 1.6 Transitions® XTRActive™	1.6	-6.00 to +6.00	6, 7, 8	Gray	X6 SV
Superlite 1.6 Polarized	1.6	-6.00 to +6.00	6, 7, 8	Gray & Brown	PZ6 SV

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.