

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Color	Invisible Marking
Hard Resin	1.498	-14.00 to +8.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8, 10	Clear	15 S FF SV
Transitions®	1.5	-14.00 to +8.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8, 10	Gray & Brown	TR S FF SV
Transitions® XTRActive™	1.5	-14.00 to +8.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8, 10	Gray	X5 S FF SV
Polarized	1.498	-14.00 to +8.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8, 10	Gray & Brown	PZ S FF SV
Drivewear®	1.5	-14.00 to +8.00	0.5, 2, 3.75, 4.25, 5.25, 6, 7, 8, 10	Drivewear®	DW S FF SV
DLC™ (Trivex)	1.53	-14.00 to +7.50	0.5, 2, 3, 4, 5, 6, 7, 8, 9	Clear	TX S FF SV
DLC™ (Trivex) Transitions®	1.53	-14.00 to +7.50	0.5, 2, 3, 4, 5, 6, 7, 8, 9	Gray & Brown	TTX S FF SV
DLC™ (Trivex) Polarized	1.53	-14.00 to +7.50	0.5, 2, 3, 4, 5, 6, 7, 8, 9	Gray & Brown	PTX S FF SV
1.56	1.56	-14.00 to +8.50	0.5, 2, 3, 4, 5, 6, 7, 8, 9	Clear	56 S FF SV
Polycarbonate	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5.25, 6.25, 7.5, 8.5	Clear	PC S FF SV
Polycarbonate Transitions®	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5.25, 6.25, 7.5, 8.5	Gray & Brown	PT S FF SV

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.

Product Availability and Markings

Material	Index	Power Range	Nominal Base Curve*	Color	Invisible Marking
Polycarbonate Transitions® XTRActive™	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5.25, 6.25, 7.5, 8.5	Gray	PX FF SV
Polycarbonate Polarized	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5.25, 6.25, 7.5, 8.5	Gray & Brown	PPZ FF SV
Polycarbonate Drivewear®	1.586	-14.00 to +8.75	0.5, 2, 3, 4, 5.25, 6.25, 7.5, 8.5	Drivewear®	PD FF SV
Superlite™ 1.6	1.6	-14.00 to +8.00	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	Clear	16 FF SV
Superlite 1.6 Transitions®	1.6	-14.00 to +8.00	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	Gray & Brown	T6 FF SV
Superlite 1.6 Transitions® XTRActive™	1.6	-14.00 to +8.00	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	Gray	X6 FF SV
Superlite 1.6 Polarized	1.6	-14.00 to +8.00	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75	Gray & Brown	PZ6 FF SV
Superlite 1.67	1.67	-16.00 to +16.00	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	Clear	67 FF SV
Superlite 1.67 Transitions®	1.67	-16.00 to +16.00	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	Gray & Brown	T67 FF SV
Superlite 1.67 Transitions® XTRActive™	1.67	-16.00 to +16.00	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	Gray	X67 FF SV
Superlite 1.67 Polarized	1.67	-16.00 to +16.00	0.5, 2, 3, 4, 4.75, 5.75, 6.75, 7.75, 9.5, 11, 12.5	Gray & Brown	PZ67 FF SV
Superlite 1.74	1.74	-20.00 to +15.50	0.5, 1, 1.5, 2.5, 3.5, 4.75, 6, 7.25, 8.25, 9, 9.75, 10.5, 11.25, 12	Clear	74 FF SV

* The true base curve can be selected from a range of curves that has been defined in the Freeform Software.

