

Shamir Relax™ The Life Improvement Lens

Lens Characteristics:

- Freeform Single Vision lens with 0.65 [D] ADD power in the lowest part of the lens
- Marked with both visible ink and invisible markings
- Aspheric design with no distortions
- Minimum fitting height: 16mm
- Fitted as regular SV lens

Uses:

- General daily use as for any other SV lens
- Suitable for any close working environment

Advantages:

- Reduces eyestrain and fatigue associated with prolonged intermediate and near vision tasks
- Customers doesn't find it different from regular SV lenses
- No distortion or swim effects

Target Market:

- 30-45 year olds engaged in prolonged intermediate and near vision tasks



Office Freeform





Shamir Relax™



The Life Improvement Lens

Minimum Fitting Height: 16mm

Includes (enter coatings and so on)

Patient price in \$ for one lens

				
Personalized SV Lens	Clear	Hi-Index	Transitions® Brown/Gray	
Lens Type & Index	Polycarbonate 1.59	SuperLite™ 1.6	Polycarbonate 1.59	SuperLite™ 1.6
Diameter				
Prescription Range	-14.0 - +8.0	-14.0 - +8.0	-14.0 - +8.0	-14.0 - +8.0
SPH				
TOR				
Shamir Glacier Plus™				
Order Codes				

Invisible Markings	
Index/Material	
	
RLX	